

# Ventilation Mode Translator

# The R. Chatburn classification

(taken from his lecture)

## 1. Determine the Control Variable

- **Volume Control**
  - Both tidal volume and inspiratory flow are preset
    - *Example: Volume Assist/Control*
- **Pressure Control**
  - Inspiratory pressure preset or proportional to effort
    - *Examples: APRV, NAVA*
- **Time Control**
  - Only inspiratory and expiratory times preset
    - *Example: HFO*

## 2. Determine the Breath Sequence

- **All breaths are spontaneous**
  - Continuous Spontaneous Ventilation (CSV)
- **Spontaneous breaths **are** possible between mandatory breaths**
  - Intermittent Mandatory Ventilation (IMV)
    1. *set mandatory rate*
    2. *mandatory breaths suppressed by spont rate*
    3. *mandatory breaths suppressed by spont MV*
- **Spontaneous breaths **are not** possible between mandatory breaths**
  - Continuous Mandatory Ventilation (CMV)

# The R. Chatburn classification

(taken from his lecture)


## 3. Determine the Target Scheme(s)

1. **Set-point** – all parameters are operator pre-set
2. **Dual** – ventilator switches between VC and PC
3. **Bio-variable** – ventilator randomly changes  $V_T$
4. **Servo** – inspiratory pressure proportional to effort
5. **Adaptive** – ventilator adjusts target with changing patient condition
6. **Optimal** – ventilator adjusts target to maximize or minimize some desired parameter
7. **Intelligent** – ventilator adjusts target using artificial intelligence tools

## Volume Assist/Control

- Operator sets inspiratory volume and flow
  - Control variable is **volume**
- Every patient effort generates a volume cycled inspiration (all breaths machine cycled)
  - Breath sequence is **CMV**
- All parameters of the control variable are operator set
  - Targeting scheme is **setpoint**

Tag is: **VC-CMV**



# Volume Controlled Modes

Chatburn Taxonomy	Hamilton Medical	Dräger	Maquet	Covidien
VC-CMV <sub>s</sub>	(S)CMV	VC-CMV VC-A/C	VC	A/C Volume
VC-CMV <sub>d</sub>	NA	VC-A/C pressure limited	VC	NA
PC-CMV <sub>a</sub>	APVcmv, CMV+	VC-CMV AutoFlow	PRVC	VC+
PC-IMV <sub>a,s</sub>	APVsimv, SIMV+	VC-SIMV Autoflow	SIMV (PRVC) + SP	SIMV-VC+
VC-IMV <sub>s,s</sub>	SIMV	VC-SIMV	SIMV(VC)+PS	SIMV
VC-IMV <sub>a,s</sub>	(is within ASV (%MinVol, Pat. height))	VC-MMV	NA	NA

# Pressure Controlled Modes

Chatburn Taxonomy	Hamilton Medical	Dräger	Maquet	Covidien
PC-CMV <sub>s</sub>	P CMV, PCV	PC-CMV PC-AC	PC	A/C Pressure
PC-IMV <sub>s,s</sub>	DuoPAP P SIMV APRV (with P <sub>support</sub> *) * Only in HAMILTON-G5	PC-BIPAP PC-SIMV PC-APRV	BI-VENT IMV(PC) + PS BI-VENT	SIMV Pressure BiLevel
PC-IMV <sub>oi,oi</sub>	ASV (%MinVol, Pat. height)	NA	NA	NA

# Pressure/Volume Support Modes

Chatburn Taxonomy	Hamilton Medical	Dräger	Maquet	Covidien
PC-CSVs	SPONT	SPN-CPAP/PS	(PS)/CPAP	SPONT
PC-CSVa	VS	SPN-CPAP/VS	VS	VS
PC-CSVr	NA	SPN-PPS	NAVA	PAV+

# Noninvasive Modes

Chatburn Taxonomy	Hamilton Medical	Dräger	Maquet	Covidien
PC-CSVs	NIV	NIV: SPN-CPAP/PS)	NIV: PS	NIV: SPONT
PC-IMVs,s	NIV-ST	NIV: PC-SIMV	NIV: SIMV(PC) + PS)	NIV: SIMV
PC-IMVs,s	nCPAP-PS <sup>1</sup>	NA	NIV: Nasal CPAP <sup>1</sup>	NIV: SPONT

<sup>1</sup>only in neonatal ventilation

# Hamilton Medical specific Modes

Chatburn Taxonomy	Hamilton Medical	Dräger	Maquet	Covidien
PC-IMVs,s	PSIMV+ (incl. Psync)	Similar: PC CMV or PC BiPAP	Similar: SIMV PC	Similar: SIMV Pressure
PC-IMVoi,oi	INTELLiVENT-ASV	NA	NA	NA

Questions to [mhimmelstoss@hamilton-medical.ch](mailto:mhimmelstoss@hamilton-medical.ch)